

dNTP Mix

Product Name	Product Code	Kit Packing
dNTP Mix, 100 mM	MBT059-1ML	1 ml
dNTP Mix, 10 mM	MBT078-1ML	1 ml

Description:

dNTP Mix, containing dATP, dCTP, dGTP, dTTP are supplied as 100 mM or 10 mM aqueous solution at pH 7.0. dNTPs are available individually as nucleotide mixes at convenient concentrations for direct use in PCR, reverse transcription, labeling and other applications.

Features:

- Greater than 99% purity confirmed by HPLC
- Free of Human and E.coli DNA
- Highly stable:
 - stable after multiple freeze-thaw cycles
 - 85-90% of dNTPs remain in triphosphate form after 30 cycles of PCR (1 min

at

94°C; 3 min at 72°C)

Applications:

For use in PCR, sequencing, nick translation, fill-in, cDNA synthesis, and TdT tailing reaction, as well as for dilution of radiolabeled dNTPs.

Storage conditions:

dNTP mix should be stored at -20°C. When stored under the recommended conditions, the product is stable for 2 years.

Quality Control:

Functionally tested in PCR
Greater than 99% purity confirmed by HPLC
Free of endo-, exo-deoxyribonuclease, ribonuclease and nicking activities

Technical Assistance

At HiMedia's, we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com

PIMBT059/78_O/0514

MBT059/78-00

Disclaimer:

User must ensure of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. Himedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research orfurther manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

🕻 🕍 A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai - 400 086, India Customer Care No.: 022-6147 1919 Email: techhelp@himedialabs.com



Registered Office :

23, Vadhani Industrial Estate,LBS Marg, Mumbai - 400 086, India.

Tel.: (022) 4017 9797 / 2500 1607 Fax: (022) 2500 2286 Commercial Office

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai - 400 086, India Tel: 00-91-22-6147 1919 Fax: 6147 1920, 2500 5764 Email: info@himedialabs.com Web: www.himedialabs.com