

DNase I Solution (1mg/ml)

Product Name DNase I Solution (1mg/ml) Product Code ML068-1ML Kit Packing 1 ml

Introduction: Deoxyribonuclease I is an RNase-free glycoprotein which specifically degrades DNA. It is a recombinant double strand specific endonuclease that catalyzes the degradation of single and double-stranded DNAs which result into di-, tri-, and oligonucleotides with 5'-phosphate and 3'-hydroxyl termini.

Description: DNase I is a glycosylated polypeptide which is commonly used for the degradation of single- and doublestranded DNAs (which are not required) into 5' phosphodinucleotide and oligonucleotide fragments. The activity of DNase I is dependent on the presence of bivalent metal ions e.g. in the presence of Mg^{2+} . Ca^{2+} or Zn^{2+} it degrades DNA by producing random single-strand nicks in the phosphate backbone and in the presence of Mn^{2+} , both strands of double-stranded DNA are cleaved leaving blunt-end fragments. Monovalent metal ions like Na⁺ and K⁺ has a negative effect on DNase I.

Application: DNase I Solution is used for specific applications where integrity of RNA has to be maintained, e.g. for the removal of genomic DNA from RNA preparations prior to RT-PCR, during nick translations, for the isolation of DNA-free RNA after *in vitro* transcription reactions, DNase etc. Moreover, it is used to catalyze random degradation of single and double-stranded DNA. DNasel solution is also used for the prevention of clumping of concentrated or cryopreserved hematopoietic cell suspensions following thawing.

Composition:

DNase I Solution (1mg/ml) consists of RNase free DNase I in a buffer which contains Sodium chloride, Magnesium chloride and Glycerol. This solution is supplied as ready to use.

Properties:

Appearance	:	Colorless solution
Clarity	:	Clear and free of particles
Bioburden	:	None detected
Suitability test	:	This solution has been tested and is suitable for use in various molecular biology applications

Storage conditions: DNase I Solution (1mg/ml) has to be stored at - 20 °C. DNase I should not be stored in frost-free freezers as temperature fluctuations will reduce its activity. Always keep DNase I on ice until ready to use. DNase I is inactivated by heating to 65°C for 10 minutes.

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at <u>mb@himedialabs.com</u>.



Registered Office :

Z3, Vadhani Industrial Estate,LBS Marg, Mumbai - 400 086, India. Tel. : (022) 4017 9797 / 2500 1607 Fax : (022) 2500 2286 Commercial Office

PIML068 0/0712

ML068-00

Tel: 00-91-22-6147 1919

Fax: 6147 1920, 2500 5764

Email : info@himedialabs.com

Web : www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein

A-516, Swastik Disha Business Park.

Via Vadhani Indl. Est., LBS Marg,

Mumbai - 400 086. India