



Technical Data

HiCrome Improved Salmonella HiVeg Agar

MV1466

HiCrome Improved Salmonella HiVeg Agar is used as an improved selective and differential medium for *Salmonella* species.

Composition**

Ingredients	Gms / Litre
HiVeg special peptone	8.000
Yeast extract	2.000
Synthetic detergent no. III	1.000
Chromogenic mixture	3.250
Agar	12.000
Final pH (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 26.25 grams in 1000 ml distilled water. Gently heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Salmonella species have been isolated from humans and almost all animals throughout the world. They cause many types of infections from mild, self-limiting gastroenteritis to life threatening typhoid fever. *Salmonella* Typhi and *Salmonella* Paratyphi A & B cause gastroenteritis, bacteremia and enteric fever, *Salmonella* Choleraesuis causes gastroenteritis and enteric fever, especially in children. *Salmonella* Typhimurium is the most frequently isolated serotype of *Salmonella* (1).

HiCrome Improved Salmonella HiVeg Agar is prepared by replacing animal based peptones with vegetable peptones and it a modification of HiCrome Improved Salmonella Agar is developed from the original formulation of Rambach (2) and is used for the differentiation of *Salmonella* species from other enteric bacteria. Rambach formulation differentiates *Salmonella* based on propylene glycol utilization and presence of a chromogenic indicator. However, HiCrome Salmonella Agar, Modified uses only a chromogenic mixture which contains chromogenic substrate and indicator dye for identification and differentiation of *Salmonella* species.

HiVeg special Peptone and yeast extract provides nitrogenous, carbonaceous compounds and other essential growth nutrients. *Escherichia coli* and *Salmonella* are easily distinguishable due to their colony characteristics. All *Salmonella* species isoalted from food or clinical sample exhibit pink to red colonies including *Salmonella* Typhi. *E. coli* exhibits a characteristic blue to purple colour, due to presence of the enzyme specific for chromogenic substrate. Synthetic detergent no. III inhibits gram-positive organisms.

Quality Control

Appearance

Cream to pinkish beige homogeneous free flowing powder

Gelling

Firm, comparable with 1.2% Agar gel.

Colour and Clarity of prepared medium

Reddish pink coloured, slightly opalescent gel forms in Petri plates

Reaction

Reaction of 2.62% w/v aqueous solution at 25°C. pH : 7.3±0.2

pH

7.10-7.50

Cultural Response

MV1466: Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Colour of Colony
<i>Bacillus subtilis</i> ATCC 6633	$\geq 10^3$	inhibited	0%	
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	$\geq 50\%$	blue to purple
<i>Proteus vulgaris</i> ATCC 13315	50-100	good	40-50%	light brown
<i>Salmonella Typhimurium</i> ATCC 14028	50-100	luxuriant	$\geq 50\%$	pink to red
<i>Salmonella Enteritidis</i> ATCC 13076	50-100	luxuriant	$\geq 50\%$	pink to red
<i>Salmonella Typhi</i> ATCC 6539	50-100	good-luxuriant	$\geq 50\%$	light pink
<i>Staphylococcus aureus</i> ATCC 25923	$\geq 10^3$	inhibited	0%	

Storage and Shelf Life

Store dehydrated powder and prepared medium at 2-8°C in tightly closed container. Use before expiry period on the label.

Reference

1. Murray P. R., Baron J. H., Pfaller M. A., Tenover J. C. and Tenover F. C., (Ed.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
2. Rambach A., 1990, Appl. Environ. Microbiol., 56:301.

Revision : 1 / 2011



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.