





Product	Code	Packing
<b>Hybridoma Tested DMSO</b>		
▽ <b>HybridoXL™ DMSO</b> Cell Culture Tested Hybridoma Tested	TC429-100ML	100ml
▽ <b>HybridoXL™ DMSO Sterile Filtered</b> Cell Culture Tested	TC433-100ML	100ml
<b>Freezing Medium for Hybridoma</b>		
* <b>CryoXL™ Hybridoma Cell Freezing Medium, DMSO, 1X</b> w/ FBS and DMSO w/o Antibiotics	TCL123-50ML	50ml
<b>Hybridoma Selection Supplements</b>		
* <b>HAT Medium Supplement 50X, Liquid</b> w/ 680.5 mg/litre Hypoxanthine, 8.8 mg/litre Aminopterin and 193.8 mg/litre Thymidine in Phosphate Buffered Saline 10ml HAT supplement is sufficient for 500ml medium	TCL072-10ML TCL072-100ML	5x10ml 10x10ml 100ml
* <b>HT Medium Supplement 50X, Liquid</b> w/ 680.5 mg/litre Hypoxanthine and 193.8 mg/litre Thymidine in Phosphate Buffered Saline 10ml HT supplement is sufficient for 500ml medium	TCL073-10ML TCL073-100ML	5x10ml 10x10ml 100ml

Product	Code	Packing
<b>Transfection Platform</b>		
<b>Media and Kits</b>		
<b>HiTransfect™ Transfection Medium, Serum Free, Animal Component Free, Low protein</b> w/ Sodium bicarbonate w/o L-Glutamine Part A - Basal Medium *Part B - Low Protein Growth Supplement	SFM015-500ML	500ml
		
<b>HiTransfect™ Transfection Medium, Reduced Serum</b> w/ Sodium bicarbonate w/o L-Glutamine Part A - Basal Medium *Part B - Reduced Serum Supplement	RSL013-500ML	500ml
		
<b>HiTransfect™ Transfection Medium, Reduced Serum</b> w/ Sodium bicarbonate w/o Phenol red and L-Glutamine Part A - Basal Medium *Part B - Reduced Serum Supplement	RSL020-500ML	500ml
		
<b>EZdetect™ β-Gal Staining Kit</b> EZdetect β Gal staining kit provides a method for determining the percentage of cells transfected with the plasmid expressing the β-galactosidase gene, LacZ. X-Gal acts as a substrate for β-Galactosidase and is hydrolyzed to form a blue coloured product, staining the transfected cells blue in colour.	CCK064-1KT	Kit for 300 tests in a 35mm dish
		

ACC