# **PTS** Pallet Truck Scale





Designed for commercial and industrial use, Adam Equipment's PTS pallet truck scale offers the durability and convenience needed in commercial and industrial settings. PTS facilitates the transport of items from one place to another in warehouses, shipping departments or on production floors. This scale is equipped with steel pallet forks, non-marring polyurethane wheels, and a zero-radius pivoting handle. The indicator is IP66 rated for protection against water and dust. With a bright LCD that is easily visible in different lighting levels, the PTS offers a wide range of features and functions to provide fast, accurate results.

# Expedite warehouse operations and reduce transportation overcharges with the PTS



#### **Features**

- Vivid, backlit LCD easily visible in any lighting conditions
- Programmable backlight can be set to "always on," "always off" or "light only when weighing"
- Capacity tracker built into display for easily monitoring possible overloads
- Checkweighing LEDs clearly show under, over and acceptable limits
- Solid non-marring polyurethane wheels
- Colour-coded keys facilitate quick recognition of the most frequently used buttons
- RS-232 interface is available to provide speedy connection to computers and printers
- Splash-proof housing, keypad and display are easy to clean
- · Wheels swivel, allowing easy movement
- Audible overload alarm warns when capacity is exceeded
- Checkweighing with audible alarm
- · External calibration allows for verification and adjustment with weights
- Parts counting with preset sample sizes
- · Auto-clear function enables the last item to be reset each time the scale returns to zero
- Selectable digital filtering for animal/dynamic weighing enables consistent results for moving subjects
- Zero-tracking feature ensures display returns to zero reading
- Rechargeable battery included for operation almost anywhere
- Programmable auto power-off to save energy
- AC adapter included

#### **Applications**

- Weighing
- Parts counting
- Percentage weighing
- Checkweighing
- Dynamic / animal weighing

#### About ADAM:

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

## **PTS** Pallet Truck Scale

 		®
	<b>AD</b> A	
	1 <i>DF</i>	

Model	PTS 2000	
Capacity	2000kg	
Readability	0.5kg	
Repeatability (S.D.)	1kg	
Linearity (+/-)	1kg	
Fork	160x1150mm	
Weighing Units	kg, lb, lb:oz, oz, N	
Stabilisation Time (sec)	2	
Interface	RS-232	
Calibration	External calibration	
Display	Backlit LCD with 40mm -high digits	
Power Supply	12VDC 50/60Hz 800mA adapter / Internal rechargeable battery	
Operating Temperature	-10° to 40°C	
Housing	Mild steel base	
Overall Dim.	550x1525x1175mm (wxdxh)	
Net Weight	125kg	

#### Accessories

Item number	Description
3102311619	In-use wet cover
700400103	. RS-232 cable
3074010507	. RS-232 to USB adapter (must order corresponding RS-232 cable)
1120011156	. ATP thermal printer
3126011263	. ATP thermal printer paper
3126011281	. ATP thermal printer paper (pack of 10)
600002028	. AdamDU data collection program
700660289	. Calibration certificate

#### **Key Accessories**



### ATP thermal printer

Easy-to-use features combined with quick set-up equals a truly versatile thermal printer. A single push-button feed key and three easy-to-read LED indicators make operation simple.



#### In-use wet cover

In-use covers provide an easy and affordable way to help keep the balance clean and protected.

### AdamDU data collection program

Adam DU software connects a balance to the laboratory network, helping to achieve seamless data collection, transfer and analysis. DU features an easy-to-use interface with a main command window, allowing users to quickly collect information from up to eight Adam balances simultaneously using RS-232 or USB connections.

Supplied by

### www.adamequipment.com