


MASS SPECTROMETRY REPORT

13-Mar-2017
18:14:31

Probe: ESI Capillary:3.00KV
Cone: 50v Extractor: 5v
Desovation Temp: 350 Gas Flow: 350
Scan ES+
8.05e5

170313-62 302 (5.582) Cm (299:303)



Sample Description

Analyzed date: 2014-03-13
Analyst: YU
Sample: C-12
M.W.: 1214.41
Lot. No.: P140303-CJ272352

Instrument

Waters ZQ2000
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN