


MASS SPECTROMETRY REPORT

14-Mar-2017
10:56:11

Probe: ESI Capillary:3.00KV
Cone: 50v Extractor: 5v
Desovation Temp: 350 Gas Flow: 350
Scan ES+
5.24e5

170314-19 247 (4.565) Cm (246:247)



Sample Description

Analyzed date: 2014-03-14
Analyst: YU
Sample: F-14
M.W.: 1742.80
Lot. No.: P140303-CJ566047

Instrument

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN