


MASS SPECTROMETRY REPORT

29-Dec-2017
11:32:49

Probe: ESI
Cone: 50v
Desolvation Temp :350

Capillary:3.00KV
Extractor: 5v
Gas Flow : 350
Scan ES+
3.53e5

171229-30 88 (1.626) Cm (88:91)



Sample Description

Analyzed date: 2016-12-29
Analyst: YU
Sample: R-16
M.W.: 1901.99
Lot. No.: P161227-WD627155

Instrument

Waters ZQ2000
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN