

HPLC REPORT

Sample: L-16

Analyzed date: 2016-12-27

Analyst: LY

Lot. No.: P161216-HX627238

Column: Vydac 218TP C18, 4.6*250mm, 5µm

Solvent A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	21%	79%
25.0min	46%	54%
25.1min	100%	0%
30.0min	Stop	

Volume: 10µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Peak No	Ret Time	Height	Area	Conc.
1	10.058	10564.694	269406.156	5.5604
2	10.220	38455.137	157700.766	3.2549
3	10.358	344879.844	4309046.500	88.9362
4	10.358	12995.277	108945.398	2.2486

Total				100.00
-------	--	--	--	--------