


HPLC REPORT

Sample: S-15 Analyzed date: 2016-1-3
Analyst: LY
Lot. No.: P161216-HX627241
Column: Kromasil-C18, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:

	A	B
0.0min	16%	84%
25.0min	41%	59%
25.1min	100%	0%
30.0min	Stop	

Volume: 10µl
Wavelength: 220nm
Flow rate: 1.0ml/min



Peak No	Ret Time	Height	Area	Conc.
1	7.987	1196.061	9828.300	0.2052
2	10.732	1136.257	17651.307	0.3686
3	12.403	962.209	4600.925	0.0961
4	12.647	251527.250	4443326.500	92.7854
5	13.832	7151.617	313415.281	6.5447
Total				100.00