


MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2016/1/3

Analyst: Shen

Sample: S-15

M.W.: 1783.91

Lot. No.: P161216-HX627241

Instrument

SHIMADZU LCMS-2020

Probe:

ESI

Probe Bias:

+ 4.5kv

Nebulizer Gas Flow:

1.5L/min

Detector:

1.2kv

CDL:

-20.0v

T. Flow:

0.2ml/min

CDL Temp.:

250 °C

B. Conc.:

50%H2O/50%ACN

Block Temp.:

400 °C