


HPLC REPORT

Sample: E-16 Analyzed date: 2016-12-28
 Analyst: LY
 Lot. No.: P161216-HX627437
 Column: Vydac 218TPC18, 4.6*250mm, 5µm
 Solvent A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	15%	85%
25.0min	40%	60%
25.1min	100%	0%
30.0min	Stop	

Volume: 10µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Peak No	Ret Time	Height	Area	Conc.
1	9.082	2691.337	25968.447	1.6913
2	10.930	381.777	3029.453	0.1973
3	12.193	1700.775	30539.844	1.9890
4	12.848	284.000	2166.905	0.1411
5	13.700	2366.078	6278.461	0.4089
6	13.892	154473.891	1446701.500	94.2227
7	13.892	3112.704	20722.139	1.3496
Total				100.00