

HPLC REPORT

Sample: V-16
Analyst: WGJ
Lot. No.: P161219-SW627439
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water


Analyzed date: 2016-12-21

Gradient:

	A	B
0.0min	11%	89%
25.0min	36%	64%
25.1min	100%	0%
30.0min		

Stop

Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	9.647	41676	1.806
2	9.803	2238407	96.98
3	10.027	27961	1.212

Total 100