

HPLC REPORT

Sample: G-16

Analyzed date: 2016-12-23

Analyst: FS

Lot. No.: P161219-SW627446

Column: 4.6*250mm, Welch CX-C18 5um

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	15%	85%
25.0min	40%	60%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Conc.	Area
1	9.352	1.25	85680
2	10.039	5.396	369789
3	10.200	87.4	5989221
4	10.424	5.958	408341
Total		100	6853031