


# HPLC REPORT

Sample: T-15 Analyzed date: 2016-12-25  
 Analyst: LLR  
 Lot. No.: P161219-SW627263  
 Column: Kromasil C18, 4.6\*250mm, 5µm  
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile  
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water  
 Gradient:
 

	A	B
0.0min	24%	76%
25.0min	49%	51%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl  
 Wavelength: 220nm  
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	11.162	6660.585	0.3072
2	11.743	14743.087	0.6800
3	12.923	2087381.125	96.2802
4	12.923	59242.277	2.7325
Total		2168027.074	100