


HPLC REPORT

Sample: G-16
Analyst: WGJ
Lot. No.: P161219-SW627462
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:
0.0min 19% 81%
25.0min 44% 56%
25.1min 100% 0%
30.0min Stop
Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min

Analyzed date: 2016-12-22



Rank	Time	Area	Conc.
1	12.598	420866.594	4.9411
2	12.882	7828125.500	91.9037
3	12.882	268756.469	3.1553
Total		100	