

HPLC REPORT

Sample: S-16

Analyzed date: 2016-12-24

Analyst: TJF

Lot. No.: P161219-CQ627468

Column: SinoChrom ODS-BP, 4.6*250mm, 5µm

Solvent A A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B B: 0.1% Trifluoroacetic Acid in 100% Water


Gradient:

	A	B
0.0min	6%	94%
25.0min	31%	69%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Conc	Area	Height
1	12.123	1.272	149697	11237
2	12.413	96.88	11405026	581516
3	12.946	1.85	217787	12203
Total		100	11772510	604956