


MASS SPECTROMETRY REPORT



Sample Description
Analyzed date: 2016/12/24
Analyst: Shen
Sample: S-16
M.W.: 1820.86
Lot. No.: P161219-CQ627468

Instrument	SHIMADZU LCMS-2020		
Probe:	ESI	Probe Bias:	+ 4.5kv
Nebulizer Gas Flow:	1.5L/min	Detector:	1.2kv
CDL:	-20.0v	T. Flow:	0.2ml/min
CDL Temp.:	250 °C	B. Conc.:	50%H2O/50%ACN
Block Temp.:	400 °C		