

HPLC REPORT

Sample: H-16P Analyzed date: 2016-12-28

Analyst: YW

Lot. No.: P161226-SY627149

Column: 4.6×250mm, Kromasil 100-5-C18

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	15%	85%
25.0min	40%	60%
25.1min	100%	0%
30.0min	Stop	

Volume: 5µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Peak	Time	Conc.	Area	Height
1	9.351	0.1863	2948	355
2	12.173	3.363	53227	4670
3	12.821	93.85	1485263	138250
4	13.175	2.302	36425	4302
5	14.302	0.2927	4632	610
Total		100	1582495	148187