


MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2016-12-28
Analyst: SHAO
Sample: H-16P
M.W.: 1905.03
Lot. No.: P161226-SY627149

Instrument

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

SHIMADZU LCMS-2020

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN