


HPLC REPORT

Sample: S-16
 Analyst: LLR
 Lot. No.: P161217-MX627387
 Column: Kromasil C18, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	8%	92%
25.0min	33%	67%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min

Analyzed date: 2016-12-23



Rank	Time	Area	Conc.
1	11.575	24004.326	1.5591
2	13.920	24118.430	1.5665
3	14.110	1459623.625	94.8055
4	14.110	31852.352	2.0689
Total		1539598.732	100