

HPLC REPORT

Sample: M-16

Analyzed date: 2016-12-24

Analyst: WGJ

Lot. No.: P161217-MX627391

Column: Kromasil C18, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	28%	72%
25.0min	53%	47%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 214nm

Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	7.550	33140	1.24
2	8.382	2235	0.08
3	8.571	2371969	88.87
4	8.964	261761	9.81
Total			100