


HPLC REPORT

Sample: Q-16
Analyst: LLR
Lot. No.: P161219-YS627071
Column: Kromasil C18, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:

Analyzed date: 2016-12-25

	A	B
0.0min	8%	92%
25.0min	33%	67%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	18.578	218172	3.894
2	19.002	5076502	90.61
3	19.350	307754	5.493
Total		5602428	100