


HPLC REPORT

Sample: L-15 Analyzed date: 2016-1-3
 Analyst: LY
 Lot. No.: P161219-YS627080
 Column: Kromasil-C18, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	17%	83%
25.0min	42%	58%
25.1min	100%	0%
30.0min	Stop	

Volume: 10µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Peak No	Ret Time	Height	Area	Conc.
1	13.802	12708.655	210559.688	6.7856
2	13.957	197258.313	2796731.750	90.1294
3	13.957	11413.982	95727.016	3.0850
Total				100.00