


HPLC REPORT

Sample: L-16 Analyzed date: 2016-12-27
Analyst: LY
Lot. No.: P161219-YS627082
Column: Vydac 218TP C18, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:

	A	B
0.0min	17%	83%
25.0min	42%	58%
25.1min	100%	0%
30.0min	Stop	

Volume: 10µl
Wavelength: 220nm
Flow rate: 1.0ml/min



Peak No	Ret Time	Height	Area	Conc.
1	9.665	743.793	4883.295	0.0589
2	13.312	6630.495	62377.344	0.7527
3	13.458	673196.688	8074000.500	97.4335
4	13.458	22927.080	136701.750	1.6497
5	15.607	1290.176	8713.293	0.1051
Total				100.00