


MASS SPECTROMETRY REPORT



Sample Description
 Analyzed date: 2016/12/26
 Analyst: Shen
 Sample: L-16
 M.W.: 2080.28
 Lot. No.: P161219-YS627082

Instrument	SHIMADZU LCMS-2020		
Probe:	ESI	Probe Bias:	+ 4.5kv
Nebulizer Gas Flow:	1.5L/min	Detector:	1.2kv
CDL:	-20.0v	T. Flow:	0.2ml/min
CDL Temp.:	250 °C	B. Conc.:	50%H2O/50%ACN
Block Temp.:	400 °C		