

HPLC REPORT

Sample: I-16

Analyzed date: 2016-12-28

Analyst: ZK

Lot. No.: P161221-CJ627101

Column: Kromasil 100-5C18, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water


Gradient :

	A	B
0.0min	10%	90%
25.0min	35%	65%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Quantity	Area	Height
1	5.508	1.4820	62358	4799
2	6.621	3.0302	127503	7720
3	6.996	88.0954	3706802	531215
4	7.260	2.1518	90543	8068
5	7.615	5.2406	220508	25979
Total		100	4207714	577781