

HPLC REPORT

Sample: H-16

Analyzed date: 2016-12-22

Analyst: ZK

Lot. No.: P161219-SY627149

Column: Kromasil 100-5C18, 4.6*250mm, 5µm

Solvent A A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient :

	A	B
0.0min	20%	80%
25.0min	45%	55%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min

mV



Rank	Time	Quantity	Area	Height
1	9.727	0.6512	31894	6964
2	9.986	4.3250	211830	17591
3	10.190	87.1730	4269530	703695
4	10.600	7.8508	384513	40959
Total		100	4897767	769209