


MASS SPECTROMETRY REPORT

22-Dec-2017
10:12:50

Probe: ESI
Cone: 50v
Desolvation Temp :350

Capillary:3.00KV
Extractor: 5v
Gas Flow : 350
Scan ES+
1.41e6

171222-11 450 (8.317)



Sample Description

Analyzed date: 2016-12-22
Analyst: YU
Sample: H-16
M.W.: 1825.05
Lot. No.: P161219-SY627149

Instrument

Waters ZQ2000
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN