


HPLC REPORT

Sample: I-16 Analyzed date: 2016-12-26
 Analyst: LY
 Lot. No.: P161215-LL627101
 Column: Vydac 218TPC18, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	16%	84%
25.0min	41%	59%
25.1min	100%	0%
30.0min	Stop	

Volume: 10µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Peak No	Ret Time	Height	Area	Conc.
1	10.890	396.482	5222.988	0.1831
2	11.268	2878.042	27514.020	0.9643
3	11.517	9726.922	59281.449	2.0777
4	11.697	22237.609	2654132.500	93.0236
5	11.697	4843.090	107031.750	3.7513
Total				100.00