


MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2016/12/26

Analyst: Shen

Sample: I-16

M.W.: 1679.80

Lot. No.: P161215-LL627101

Instrument

SHIMADZU LCMS-2020

Probe:

ESI

Probe Bias:

+ 4.5kv

Nebulizer Gas Flow: 1.5L/min

Detector:

1.2kv

CDL: -20.0v

T. Flow:

0.2ml/min

CDL Temp.: 250 °C

B. Conc.:

50%H2O/50%ACN

Block Temp.: 400 °C