


MASS SPECTROMETRY REPORT

6.82e7



Sample Description

Analyzed date: 2013-07-16
 Analyst: YU
 Sample: 5-333 CE-15
 M.W.: 1719.81
 Lot. No.: P130704-LW327169

Instrument

Waters ZQ2000
 Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp.: 250°C
 Block Temp.: 200°C

Probe Bias: + 4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. Conc.: 50%H2O/50% ACN