

HPLC REPORT

Sample: 5-302 SD-15

Analyzed date: 2013-07-24

Analyst: ZXX

Lot. No.: P130705-XZ358236

Column: Kromasil C18, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	10%	90%
25.0min	35%	65%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	12.103	88301.281	2.2212
2	12.622	3769712.250	94.8270
3	12.622	105133.375	2.6446
4	15.847	12211.601	0.3072
Total		3975358.507	100