

HPLC REPORT

Sample: 5-345 VV-18

Analyzed date: 2013-07-17

Analyst: LYZ

Lot. No.: P130704-LW272428

Column: SinoChrom ODS-BP, 4.6*250mm, 5 μ m

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	24%	76%
25.0min	49%	51%
25.1min	100%	0%
30.0min	Stop	

Volume: 20 μ l

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	10.232	11016.130	0.5705
2	10.505	1718146.000	88.9722
3	10.505	201943.203	10.4574
Total		1931105.333	100