

HPLC REPORT

Sample: 5-354 CL-20

Analyzed date: 2013-07-19

Analyst: LYZ

Lot. No.: P130704-LW358278

Column: SinoChrom ODS-BP, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	20%	80%
25.0min	45%	55%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	4.362	40219.246	1.7736
2	4.977	10055.001	0.4434
3	5.948	6716.596	0.2962
4	8.052	15841.994	0.6986
5	11.248	30415.123	1.3413
6	11.463	1960259.875	86.4439
7	11.852	41758.230	1.8415
8	12.287	162402.359	7.1616
Total		2267668.425	100