


MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2013-7-19
Analyst: FENG
Sample: 5-354 CL-20
M.W.: 2562.95
Lot. No.: P130704-LW358278

Instrument

SHIMADZU LCMS-2010EV
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C
Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN