


HPLC REPORT

Sample: 5-388 CQ-16 Analyzed date: 2013-07-19
 Analyst: LYZ
 Lot. No.: P130704-LW272430
 Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	20%	80%
25.0min	45%	55%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	4.857	12417.501	0.5879
2	10.917	6151.728	0.2913
3	11.073	1980991.500	93.7930
4	11.073	73947.430	3.5012
5	11.642	38581.488	1.8267
Total		2112089.647	100