

HPLC REPORT

Sample: 5-398 CL-15

Analyzed date: 2013-07-23

Analyst: LYZ

Lot. No.: P130704-LW358310

Column: SinoChrom ODS-BP, 4.6*250mm, 5 μ m

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	15%	85%
25.0min	40%	60%
25.1min	100%	0%
30.0min	Stop	

Volume: 20 μ l

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	9.780	68916.883	3.0158
2	10.250	63544.203	2.7807
3	10.472	2049138.750	89.6696
4	10.472	91542.078	4.0059
5	11.258	12067.371	0.5281

Total 2285209.285 100