

HPLC REPORT

Sample: 5-402 CR-15

Analyzed date: 2013-07-23

Analyst: LYZ

Lot. No.: P130704-LW358314

Column: SinoChrom ODS-BP, 4.6*250mm, 5 μ m

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	20%	80%
25.0min	45%	55%
25.1min	100%	0%
30.0min	Stop	

Volume: 20 μ l

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	4.617	11534.751	1.3517
2	5.810	4572.400	0.5358
3	7.783	50632.008	5.9331
4	9.690	5574.832	0.6533
5	9.913	761109.000	89.1878
6	10.198	19955.256	2.3384
Total		853378.247	100