

HPLC REPORT

Sample: 5-409 CD-15

Analyzed date: 2013-07-25

Analyst: LYZ

Lot. No.: P130704-LW272445

Column: SinoChrom ODS-BP, 4.6*250mm, 5 μ m

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	14%	86%
25.0min	39%	61%
25.1min	100%	0%
30.0min	Stop	

Volume: 20 μ l

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	13.378	40764.371	1.7704
2	13.540	2222342.500	96.5154
3	13.540	33017.027	1.4339
4	15.718	6453.555	0.2803
Total		2302577.453	100