


# HPLC REPORT

Sample: 5-404 CT-17 Analyzed date: 2013-07-23  
Analyst: LYZ  
Lot. No.: P130704-LW358316  
Column: SinoChrom ODS-BP, 4.6\*250mm, 5µm  
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile  
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water  
Gradient:  
  A  B  
          0.0min          12%  88%  
          25.0min         37%  63%  
          25.1min         100%   0%  
          30.0min  Stop  
Volume: 20µl  
Wavelength: 220nm  
Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	8.158	7589.801	0.4604
2	10.637	17332.395	1.0514
3	10.793	1518789.750	92.1295
4	11.083	104826.617	6.3588
Total		1648538.536	100