


# MASS SPECTROMETRY REPORT

1.84e7



Sample Description

Analyzed date: 2013-07-11  
 Analyst: YU  
 Sample: 5-297 CD-8  
 M.W.: 998.14  
 Lot. No.: P130704-ML358233

Instrument

Waters ZQ2000  
 Probe: ESI  
 Nebulizer Gas Flow: 1.5L/min  
 CDL: -20.0v  
 CDL Temp.: 250°C  
 Block Temp.: 200°C

Waters ZQ2000

Probe Bias: + 4.5kv  
 Detector: 1.5kv  
 T. Flow: 0.2ml/min  
 B. Conc.: 50% H2O/50% ACN