


MASS SPECTROMETRY REPORT



1.34e8

Sample Description

Analyzed date: 2013-07-10
 Analyst: YU
 Sample: 6-018N CE-10
 M.W.: 1234.31
 Lot. No.: P130704-LW249021

Instrument

Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp.: 250 °C
 Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: + 4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. Conc.: 50%H2O/50%ACN