

HPLC REPORT

Sample: 6-018P CE-10

Analyzed date: 2013-07-11

Analyst: LYZ

Lot. No.: P130704-LW329009

Column: SinoChrom ODS-BP, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water


Gradient:

	A	B
0.0min	10%	90%
25.0min	35%	65%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	4.553	50815.797	1.3557
2	11.773	3436848.000	91.6921
3	11.773	260586.219	6.9522
Total		3748250.016	100