


HPLC REPORT

Sample: C-12P
Analyst: YW
Lot. No.: P140929-SY608867
Column: 4.6×250mm, Kromasil 100-5-C18
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:

	A	B
0.0min	15%	85%
25.0min	40%	60%
25.1min	100%	0%
30.0min	Stop	

Volume: 5µl
Wavelength: 220nm
Flow rate: 1.0ml/min

Analyzed date: 2014-10-13



Peak	Time	Height	Area	Conc.
1	9.378	10633.465	76550.539	1.7903
2	9.638	349403.938	3803865.750	88.9592
3	10.012	16552.594	188002.609	4.3967
4	10.318	5944.717	63008.258	1.4735
5	10.812	3167.134	68712.500	1.6069
6	18.940	433.607	3593.601	0.0840
7	22.087	510.069	6054.454	0.1416
8	22.387	923.013	7516.280	0.1758
9	22.687	1304.219	14495.018	0.3390
10	22.942	2024.866	17127.996	0.4006
11	24.150	2021.926	27040.000	0.6324
Total				100.000