


MASS SPECTROMETRY REPORT

11-Oct-2017
10:48:28

Probe: ESI
Cone: 50v
Desovation Temp :350
Capillary:3.00KV
Extractor: 5v
Gas Flow : 350
Scan ES-
1.87e4

171011-23- 122 (2.255) Cm (122:127)



Sample Description

Analyzed date: 2014-10-11
Analyst: YU
Sample: C-16P
M.W.: 1431.41
Lot. No.: P140921-HS609335

Instrument

Waters ZQ2000
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN