

HPLC REPORT

Sample: K-9

Analyzed date: 2014-10-19

Analyst: ZK

Lot. No.: P141010-WD612010

Column: Kromasil 100-5C18, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient :


	A	B
0.0min	20%	80%
25.0min	45%	55%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min

mV



Rank	Time	Quantity	Area	Height
1	10.879	0.6965	33841	3303
2	11.023	95.9004	4659498	543396
3	11.311	2.0658	100370	11262
4	11.922	1.3373	64974	8315
Total		100	4858683	566276