


HPLC REPORT

Sample: T-14
 Analyst: WGJ
 Lot. No.: P160922--YJ608754
 Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	13%	87%
25.0min	38%	62%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min

Analyzed date: 2016-11-29



Rank	Time	Area	Conc.
1	2.767	170414.219	2.2714
2	3.377	4338.308	0.0578
3	3.550	4412315.000	58.8110
4	3.942	646375.438	8.6154
5	4.518	13552.000	0.1806
6	5.085	20971.752	0.2795
7	7.485	62359.805	0.8312
8	8.042	902424.500	12.0283
9	9.457	51330.012	0.6842
10	12.115	1218457.375	16.2406

Total 100