


MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2014/10/9

Analyst: Shen

Sample: C-14P

M.W.: 1688.74

Lot. No.: P140922--YJ608771

Instrument

SHIMADZU LCMS-2020

Probe:

ESI

Probe Bias:

+ 4.5kv

Nebulizer Gas Flow: 1.5L/min

CDL: -20.0v

CDL Temp.: 250 °C

Block Temp.: 400 °C

Detector: 1.2kv

T. Flow: 0.2ml/min

B. Conc.: 50%H2O/50%ACN