

HPLC REPORT


Sample: C-14P
Analyst: WGJ
Lot. No.: P140922--YJ608772
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Analyzed date: 2014-09-25

Gradient:

	A	B
0.0min	17%	83%
25.0min	42%	58%
25.1min	100%	0%
30.0min		Stop

Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min
mV



Rank	Time	Area	Conc.
1	12.050	55753	1.543
2	12.156	3466262	95.95
3	12.500	90648	2.509
Total			100